

ZSCAN26 Antibody



PACO44734

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZSCAN26

Uniprot

Q16670

Synonyms:

Zinc finger and SCAN domain-containing protein 26 (Protein SRE-ZBP) (Zinc finger protein 187), ZSCAN26, ZNF187


Immunogen:

Recombinant Human Zinc finger and SCAN domain-containing protein 26 protein (239-478AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO44734 at dilution of 1:100.