

ZNF20 Antibody



PACO44737

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF20

Uniprot

P17024

Synonyms:

Zinc finger protein 20 (Zinc finger protein KOX13), ZNF20, KOX13


Immunogen:

Recombinant Human Zinc finger protein 20 protein (1-200AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: ZNF20 antibody at 2.53 μ g/ml + MCF-7 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 62 kDa. Observed band size: 62, 40 kDa.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44737 at dilution of 1:100.