

ZNF264 Antibody



PACO44741

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF264

Uniprot

O43296

Synonyms:

Zinc finger protein 264, ZNF264, KIAA0412

Immunogen:

Recombinant Human Zinc finger protein 264 protein (1-220AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: ZNF264 antibody at 1.02 μ g/ml + Jurkat whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 71 kDa. Observed band size: 71 kDa.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44741 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO44741 at dilution of 1:100.