

ZNF544 Antibody



PACO44744

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF544

Uniprot

Q6NX49

Synonyms:

Zinc finger protein 544, ZNF544


Immunogen:

Recombinant Human Zinc finger protein 544 protein (1-230AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO44744 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44744 at dilution of 1:100.