

CXXC4 Antibody



PACO44757

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Acts as a negative regulator of the Wnt signaling pathway via its interaction with DVL1.

Gene ID:

CXXC4

Uniprot

Q9H2H0

Synonyms:

CXXC-type zinc finger protein 4 (Inhibition of the Dvl and axin complex protein), CXXC4, IDAX


Immunogen:

Recombinant Human CXXC-type zinc finger protein 4 protein (1-198AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO44757 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44757 at dilution of 1:100.