

PACO44772

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytosol, inositol-1,4,5-trisphosphate 3-kinase activity, cellular response to calcium ion, inositol phosphate metabolic process, inositol trisphosphate metabolic process, signal transduction.

Gene ID:

ITPKB

Uniprot

P27987

Synonyms:

Inositol-trisphosphate 3-kinase B (EC 2.7.1.127) (Inositol 1,4,5-trisphosphate 3-kinase B) (IP3 3-kinase B) (IP3K B) (InsP 3-kinase B), ITPKB

Immunogen:

Recombinant Human Inositol-trisphosphate 3-kinase B protein (1-260AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A