

TNNI3K Antibody



PACO44805

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May play a role in cardiac physiology.

Gene ID:

TNNI3K

Uniprot

Q59H18

Synonyms:

Serine/threonine-protein kinase TNNI3K (EC 2.7.11.1) (Cardiac ankyrin repeat kinase) (Cardiac troponin I-interacting kinase) (TNNI3-interacting kinase), TNNI3K, CARK


Immunogen:

Recombinant Human Serine/threonine-protein kinase TNNI3K protein (1-75AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO44805 at dilution of 1:100.