

# TSSK3 Antibody



PACO44807

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May be involved in a signaling pathway during male germ cell development or mature sperm function.

**Gene ID:**

TSSK3

**Uniprot**

Q96PN8

**Synonyms:**

Testis-specific serine/threonine-protein kinase 3 (TSK-3) (TSSK-3) (Testis-specific kinase 3) (EC 2.7.11.1) (Serine/threonine-protein kinase 22C), TSSK3, SPOGA3 STK22C

**Immunogen:**


Recombinant Human Testis-specific serine/threonine-protein kinase 3 protein (1-268AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human lung cancer using PACO44807 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver cancer using PACO44807 at dilution of 1:100.