

PACO44884

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Potent chemoattractant for neutrophils, and weaker for dendritic cells. Not chemotactic for T-cells, B-cells, monocytes, natural killer cells or granulocytes. Does not inhibit proliferation of myeloid progenitors in colony formation assays.

Gene ID:

CXCL14

Uniprot

O95715

Synonyms:

C-X-C motif chemokine 14 (Chemokine BRAK) (MIP-2G) (Small-inducible cytokine B14), CXCL14, MIP2G NJAC SCYB14


Immunogen:

Recombinant Human C-X-C motif chemokine 14 protein (35-111AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO44884 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO44884 at dilution of 1:100.