

CST8 Antibody



PACO44900

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Performs a specialized role during sperm development and maturation.

Gene ID:

CST8

Uniprot

O60676

Synonyms:

Cystatin-8 (Cystatin-related epididymal spermatogenic protein), CST8, CRES

Immunogen:

Recombinant Human Cystatin-8 protein (53-142AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A