

PACO44974

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Thought to be involved during neural development in axonal navigation at the ventral midline of the neural tube. In spinal chord development plays a role in guiding commissural axons probably by preventing premature sensitivity to Slit proteins thus inhibiting Slit signaling through ROBO1. Required for hindbrain axon midline crossing.

Gene ID:

ROBO3

Uniprot

Q96MS0

Synonyms:

Roundabout homolog 3 (Roundabout-like protein 3), ROBO3


Immunogen:

Recombinant Human Roundabout homolog 3 protein (1237-1386AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO44974 at dilution of 1:100.