

PACO45023

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Catalyzes the removal of sialic acid, (N-acetylneuraminic acid, moities from glycoproteins, oligosaccharides and gangliosides.

**Gene ID:**

NEU2

**Uniprot**

Q9Y3R4

**Synonyms:**

Sialidase-2 (EC 3.2.1.18) (Cytosolic sialidase) (N-acetyl-alpha-neuraminidase 2), NEU2

**Immunogen:**


Recombinant Human Sialidase-2 protein (171-380AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO45023 at dilution of 1:100.