

RNF13 Antibody



PACO45199

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

E3 ubiquitin-protein ligase that may play a role in controlling cell proliferation.

Gene ID:

RNF13

Uniprot

O43567

Synonyms:

E3 ubiquitin-protein ligase RNF13 (EC 2.3.2.27) (RING finger protein 13) (RING-type E3 ubiquitin transferase RNF13), RNF13, RZF


Immunogen:

Recombinant Human E3 ubiquitin-protein ligase RNF13 protein (204-381AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: RNF13 antibody at 1.22 μ g/ml + MCF-7 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 43, 30 kDa. Observed band size: 43 kDa.



Immunohistochemistry of paraffin-embedded human glioma using PACO45199 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO45199 at dilution of 1:100.