

RAP2B Antibody



PACO45257

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Small GTP-binding protein which cycles between a GDP-bound inactive and a GTP-bound active form. Involved in EGFR and CHRM3 signaling pathways through stimulation of PLCE1. May play a role in cytoskeletal rearrangements and regulate cell spreading through activation of the effector TNIK. May regulate membrane vesiculation in red blood cells.

Gene ID:

RAP2B

Uniprot

P61225

Synonyms:

Ras-related protein Rap-2b, RAP2B


Immunogen:

Recombinant Human Ras-related protein Rap-2b protein (1-183AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO45257 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO45257 at dilution of 1:100.