

ADAMTSL4 Antibody



PACO45279

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Positive regulation of apoptosis. May facilitate FBN1 microfibril biogenesis.

Gene ID:

ADAMTSL4

Uniprot

Q6UY14

Synonyms:

ADAMTS-like protein 4 (ADAMTSL-4) (Thrombospondin repeat-containing protein 1), ADAMTSL4, TSRC1


Immunogen:

Recombinant Human ADAMTS-like protein 4 protein (795-1074AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO45279 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO45279 at dilution of 1:100.