

# ADAM19 Antibody



PACO45407

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Participates in the proteolytic processing of beta-type neuregulin isoforms which are involved in neurogenesis and synaptogenesis, suggesting a regulatory role in glial cell. Also cleaves alpha-2 macroglobulin. May be involved in osteoblast differentiation and/or osteoblast activity in bone.

**Gene ID:**

ADAM19

**Uniprot**

Q9H013

**Synonyms:**

Disintegrin and metalloproteinase domain-containing protein 19 (ADAM 19) (EC 3.4.24.-) (Meltrin-beta) (Metalloprotease and disintegrin dendritic antigen marker) (MADDAM), ADAM19, MLTNB

**Immunogen:**


Recombinant Human Disintegrin and metalloproteinase domain-containing protein 19 protein (203-460AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO45407 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO45407 at dilution of 1:100.