

CD53 Antibody



PACO45414

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Required for efficient formation of myofibers in regenerating muscle at the level of cell fusion. May be involved in growth regulation in hematopoietic cells.

Gene ID:

CD53

Uniprot

P19397

Synonyms:

Leukocyte surface antigen CD53 (Cell surface glycoprotein CD53) (Tetraspanin-25) (Tspan-25) (CD antigen CD53), CD53, MOX44 TSPAN25


Immunogen:

Recombinant Human Leukocyte surface antigen CD53 protein (100-180AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO45414 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO45414 at dilution of 1:100.