

NRXN1 Antibody



PACO45461

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

NRXN1

Uniprot

Q9ULB1

Synonyms:

Neurexin-1 (Neurexin I-alpha) (Neurexin-1-alpha), NRXN1, KIAA0578


Immunogen:

Recombinant Human Neurexin-1 protein (31-310AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: NRXN1 antibody at 1.57 μ g/ml. Lane 1: Jurkat whole cell lysate. Lane 2: Hela whole cell lysate. Lane 3: HepG2 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 162, 165, 170, 47, 16, 51 kDa. Observed band size: 162 kDa.