

ANKS1B Antibody



PACO45508

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Isoform 2 may participate in the regulation of nucleoplasmic coilin protein interactions in neuronal and transformed cells.

Gene ID:

ANKS1B

Uniprot

Q7Z6G8

Synonyms:

Ankyrin repeat and sterile alpha motif domain-containing protein 1B (Amyloid-beta protein intracellular domain-associated protein 1) (AIDA-1) (E2A-PBX1-associated protein) (EB-1), ANKS1B


Immunogen:

Recombinant Human Ankyrin repeat and sterile α motif domain-containing protein 1B protein (999-1248AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO45508 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO45508 at dilution of 1:100.