

MRPL12 Antibody



PACO45519

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

mitochondrial inner membrane, mitochondrial large ribosomal subunit, mitochondrion, RNA binding, structural constituent of ribosome, mitochondrial translational elongation, mitochondrial translational termination, positive regulation of transcription, DNA-templated, transcription from mitochondrial promoter.

Gene ID:

MRPL12

Uniprot

P52815

Synonyms:

39S ribosomal protein L12, mitochondrial (L12mt) (MRP-L12) (5c5-2) (Mitochondrial large ribosomal subunit protein bL12m), MRPL12, MRPL7 RPML12


Immunogen:

Recombinant Human 39S ribosomal protein L12, mitochondrial protein (1-198AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO45519 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO45519 at dilution of 1:100.