

ARHGEF3 Antibody



PACO45532

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Acts as guanine nucleotide exchange factor (GEF) for RhoA and RhoB GTPases.

Gene ID:

ARHGEF3

Uniprot

Q9NR81

Synonyms:

Rho guanine nucleotide exchange factor 3 (Exchange factor found in platelets and leukemic and neuronal tissues) (XPLN), ARHGEF3


Immunogen:

Recombinant Human Rho guanine nucleotide exchange factor 3 protein (287-526AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO45532 at dilution of 1:100.