

PMPCA Antibody



PACO45552

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Cleaves presequences (transit peptides) from mitochondrial protein precursors.

Gene ID:

PMPCA

Uniprot

Q10713

Synonyms:

Mitochondrial-processing peptidase subunit alpha (EC 3.4.24.64) (Alpha-MPP) (P-55), PMPCA, INPP5E KIAA0123 MPPA


Immunogen:

Recombinant Human Mitochondrial-processing peptidase subunit α protein (456-525AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO45552 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO45552 at dilution of 1:100.