

# PTPN4 Antibody



PACO45553

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May act at junctions between the membrane and the cytoskeleton.

**Gene ID:**

PTPN4

**Uniprot**

P29074

**Synonyms:**

Tyrosine-protein phosphatase non-receptor type 4 (EC 3.1.3.48) (Protein-tyrosine phosphatase MEG1) (MEG) (PTPase-MEG1), PTPN4

**Immunogen:**


Recombinant Human Tyrosine-protein phosphatase non-receptor type 4 protein (1-240AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO45553 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO45553 at dilution of 1:100.