

PACO45646

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probably catalyzes the addition of a phosphoethanolamine moiety to lipid A. Phosphoethanolamine modification of lipid A gives polymyxin resistance.

Gene ID:

mcr1

Uniprot

A0A0R6L508

Synonyms:

Probable phosphatidylethanolamine transferase Mcr-1 (EC 2.7. -. -) (Polymyxin resistance protein MCR-1), mcr1, mcr-1

Immunogen:

Recombinant Escherichia coli Probable phosphatidylethanolamine transferase Mcr-1 protein (1-541AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A