

ADORA2A Antibody, Biotin conjugated



PACO45653

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor for adenosine. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Gene ID:

ADORA2A

Uniprot

P29274

Synonyms:

Adenosine receptor A2a, ADORA2A, ADORA2

Immunogen:

Recombinant Human Adenosine receptor A2a protein (291-412AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A