

# BUB1 Antibody, FITC conjugated



PACO45700

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

BUB1

**Uniprot**

O43683

**Synonyms:**

Mitotic checkpoint serine/threonine-protein kinase BUB1 (hBUB1) (EC 2.7.11.1) (BUB1A), BUB1, BUB1L

**Immunogen:**

Recombinant Human Mitotic checkpoint serine/threonine-protein kinase BUB1 protein (256-577AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A