

BUB1 Antibody, Biotin conjugated



PACO45701

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

BUB1

Uniprot

O43683

Synonyms:

Mitotic checkpoint serine/threonine-protein kinase BUB1 (hBUB1) (EC 2.7.11.1) (BUB1A), BUB1, BUB1L

Immunogen:

Recombinant Human Mitotic checkpoint serine/threonine-protein kinase BUB1 protein (256-577AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A