

CD83 Antibody



PACO45722

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

May play a significant role in antigen presentation or the cellular interactions that follow lymphocyte activation.

Gene ID:

CD83

Uniprot

Q01151

Synonyms:

CD83 antigen (hCD83) (B-cell activation protein) (Cell surface protein HB15) (CD antigen CD83), CD83


Immunogen:

Recombinant Human CD83 antigen protein (20-144AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO45722 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).