

EFCAB7 Antibody, FITC conjugated



PACO45796

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

primary cilium.

Gene ID:

EFCAB7

Uniprot

A8K855

Synonyms:

EF-hand calcium-binding domain-containing protein 7, EFCAB7, KIAA1799

Immunogen:

Recombinant Human EF-hand calcium-binding domain-containing protein 7 protein (1-158AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A