

# EFCAB7 Antibody, Biotin conjugated



PACO45797

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

primary cilium.

**Gene ID:**

EFCAB7

**Uniprot**

A8K855

**Synonyms:**

EF-hand calcium-binding domain-containing protein 7, EFCAB7, KIAA1799

**Immunogen:**

Recombinant Human EF-hand calcium-binding domain-containing protein 7 protein (1-158AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A