

FABP7 Antibody, HRP conjugated



PACO45823

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

B-FABP could be involved in the transport of a so far unknown hydrophobic ligand with potential morphogenic activity during CNS development. It is required for the establishment of the radial glial fiber system in developing brain, a system that is necessary for the migration of immature neurons to establish cortical layers.

Gene ID:

FABP7

Uniprot

O15540

Synonyms:

Fatty acid, binding protein, brain (Brain lipid-binding protein) (BLBP) (Brain-type fatty acid, binding protein) (B-FABP) (Fatty acid, binding protein 7) (Mammary-derived growth inhibitor related), FABP7, BLBP FABPB MRG

Immunogen:

Recombinant Human Fatty acid, binding protein, brain protein (2-132AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A