

# RIPOR2 Antibody



PACO45826

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Required for hearing. Involved in skeletal muscle development.

**Gene ID:**

RIPOR2

**Uniprot**

Q9Y4F9

**Synonyms:**

Rho family-interacting cell polarization regulator 2, RIPOR2, C6orf32 DIFF48 FAM65B KIAA0386 PL48

**Immunogen:**


Recombinant Human Rho family-interacting cell polarization regulator 2 protein (1-250AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO45826 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO45826 at dilution of 1:100.