

HFM1 Antibody, Biotin conjugated



PACO45889

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Required for crossover formation and complete synapsis of homologous chromosomes during meiosis.

Gene ID:

HFM1

Uniprot

A2PYH4

Synonyms:

Probable ATP-dependent DNA helicase HFM1 (EC 3.6.4.12) (SEC63 domain-containing protein 1), HFM1, SEC3D1

Immunogen:

Recombinant Human Probable ATP-dependent DNA helicase HFM1 protein (121-420AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A