

HS6ST1 Antibody, FITC conjugated



PACO45916

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

6-O-sulfation enzyme which catalyzes the transfer of sulfate from 3'-phosphoadenosine 5'-phosphosulfate (PAPS) to position 6 of the N-sulfoglucosamine residue (GlcNS) of heparan sulfate. Critical for normal neuronal development where it may play a role in neuron branching. May also play a role in limb development. May prefer iduronic acid, .

Gene ID:

HS6ST1

Uniprot

O60243

Synonyms:

Heparan-sulfate 6-O-sulfotransferase 1 (HS6ST-1) (EC 2.8.2. -), HS6ST1, HS6ST

Immunogen:

Recombinant Human Heparan-sulfate 6-O-sulfotransferase 1 protein (3-133AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A