

# LATS1 Antibody, HRP conjugated



PACO45971

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LATS1

**Uniprot**

O95835

**Synonyms:**

Serine/threonine-protein kinase LATS1 (EC 2.7.11.1) (Large tumor suppressor homolog 1) (WARTS protein kinase) (h-warts), LATS1, WARTS

**Immunogen:**

Recombinant Human Serine/threonine-protein kinase LATS1 protein (1-280AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A