

LMO2 Antibody, FITC conjugated



PACO45984

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Acts with TAL1/SCL to regulate red blood cell development. Also acts with LDB1 to maintain erythroid precursors in an immature state.

Gene ID:

LMO2

Uniprot

P25791

Synonyms:

Rhombotin-2 (Cysteine-rich protein TTG-2) (LIM domain only protein 2) (LMO-2) (T-cell translocation protein 2), LMO2, RBTN2 RBTNL1 RHOM2 TTG2

Immunogen:

Recombinant Human Rhombotin-2 protein (1-100AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A