

ME1 Antibody, HRP conjugated



PACO46015

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytosol, ADP binding, electron carrier activity, malate dehydrogenase (decarboxylating) (NADP+) activity, malic enzyme activity, manganese ion binding, NAD binding, NADP binding, carbohydrate metabolic process, malate metabolic process.

Gene ID:

ME1

Uniprot

P48163

Synonyms:

NADP-dependent malic enzyme (NADP-ME) (EC 1.1.1.40) (Malic enzyme 1), ME1

Immunogen:

Recombinant Human NADP-dependent malic enzyme protein (161-519AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A