

NFAT5 Antibody, Biotin conjugated



PACO46073

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcription factor involved in the transcriptional regulation of osmoprotective and inflammatory genes. Regulates hypertonicity-induced cellular accumulation of osmolytes.

Gene ID:

NFAT5

Uniprot

O94916

Synonyms:

Nuclear factor of activated T-cells 5 (NF-AT5) (T-cell transcription factor NFAT5) (Tonicity-responsive enhancer-binding protein) (TonE-binding protein) (TonEBP), NFAT5, KIAA0827 TONEBP

Immunogen:

Recombinant Human Nuclear factor of activated T-cells 5 protein (1422-1531AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A