

NR2F2 Antibody, HRP conjugated



PACO46087

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Ligand-activated transcription factor. Activated by high concentrations of 9-cis-retinoic acid, and all-trans-retinoic acid, but not by dexamethasone, cortisol or progesterone (in vitro). Regulation of the apolipoprotein A-I gene transcription. Binds to DNA site A.

Gene ID:

NR2F2

Uniprot

P24468

Synonyms:

COUP transcription factor 2 (COUP-TF2) (Apolipoprotein A-I regulatory protein 1) (ARP-1) (COUP transcription factor II) (COUP-TF II) (Nuclear receptor subfamily 2 group F member 2), NR2F2, ARP1 TFCOUP2

Immunogen:

Recombinant Human COUP transcription factor 2 protein (123-383AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A