

PCF11 Antibody



PACO46098

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Component of pre-mRNA cleavage complex II.

Gene ID:

PCF11

Uniprot

O94913

Synonyms:

Pre-mRNA cleavage complex 2 protein Pcf11 (Pre-mRNA cleavage complex II protein Pcf11), PCF11, KIAA0824


Immunogen:

Recombinant Human Pre-mRNA cleavage complex 2 protein Pcf11 protein (409-604AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO46098 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO46098 at dilution of 1:100.