

# PCF11 Antibody, HRP conjugated



PACO46099

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Component of pre-mRNA cleavage complex II.

**Gene ID:**

PCF11

**Uniprot**

O94913

**Synonyms:**

Pre-mRNA cleavage complex 2 protein Pcf11 (Pre-mRNA cleavage complex II protein Pcf11), PCF11, KIAA0824

**Immunogen:**

Recombinant Human Pre-mRNA cleavage complex 2 protein Pcf11 protein (409-604AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A