

# PHC1 Antibody, Biotin conjugated



PACO46113

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development. PcG PRC1 complex acts via chromatin remodeling and modification of histones; it mediates monoubiquitination of histone H2A 'Lys-119', rendering chromatin heritably changed in its expressibility. Required for proper control of cellular levels of GMNN expression.

**Gene ID:**

PHC1

**Uniprot**

P78364

**Synonyms:**

Polyhomeotic-like protein 1 (hPH1) (Early development regulatory protein 1), PHC1, EDR1 PH1

**Immunogen:**

Recombinant Human Polyhomeotic-like protein 1 protein (142-225AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A