

RANBP1 Antibody, Biotin conjugated



PACO46157

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Inhibits GTP exchange on Ran. Forms a Ran-GTP-RANBP1 trimeric complex. Increase GTP hydrolysis induced by the Ran GTPase activating protein RANGAP1. May act in an intracellular signaling pathway which may control the progression through the cell cycle by regulating the transport of protein and nucleic acid, across the nuclear membrane.

Gene ID:

RANBP1

Uniprot

P43487

Synonyms:

Ran-specific GTPase-activating protein (Ran-binding protein 1) (RanBP1), RANBP1

Immunogen:

Recombinant Human Ran-specific GTPase-activating protein (52-111AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A