

# SUV39H1 Antibody, HRP conjugated



PACO46259

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SUV39H1

**Uniprot**

O43463

**Synonyms:**

Histone-lysine N-methyltransferase SUV39H1 (EC 2.1.1.43) (Histone H3-K9 methyltransferase 1) (H3-K9-HMTase 1) (Lysine N-methyltransferase 1A) (Position-effect variegation 3-9 homolog) (Suppressor of variegation 3-9 homolog 1) (Su(var)3-9 homolog 1), SUV39H1, KMT1A SUV39H

**Immunogen:**

Recombinant Human Histone-lysine N-methyltransferase SUV39H1 protein (179-412AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A