

SUV39H1 Antibody, Biotin conjugated



PACO46261

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SUV39H1

Uniprot

O43463

Synonyms:

Histone-lysine N-methyltransferase SUV39H1 (EC 2.1.1.43) (Histone H3-K9 methyltransferase 1) (H3-K9-HMTase 1) (Lysine N-methyltransferase 1A) (Position-effect variegation 3-9 homolog) (Suppressor of variegation 3-9 homolog 1) (Su(var)3-9 homolog 1), SUV39H1, KMT1A SUV39H

Immunogen:

Recombinant Human Histone-lysine N-methyltransferase SUV39H1 protein (179-412AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A