

# TUBB4A Antibody, HRP conjugated



PACO46307

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

**Gene ID:**

TUBB4A

**Uniprot**

P04350

**Synonyms:**

Tubulin beta-4A chain (Tubulin 5 beta) (Tubulin beta-4 chain), TUBB4A, TUBB4 TUBB5

**Immunogen:**

Recombinant Human Tubulin beta-4A chain protein (1-444AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A