

WRB Antibody, HRP conjugated



PACO46327

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor for ASNA1/TRC40-mediated insertion of tail-anchored (TA) proteins into the ER membrane.

Gene ID:

WRB

Uniprot

O00258

Synonyms:

Tail-anchored protein insertion receptor WRB (Congenital heart disease 5 protein) (Tryptophan-rich basic protein) (WRB), WRB, CHD5

Immunogen:

Recombinant Human Tail-anchored protein insertion receptor WRB protein (38-99AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A