

# WRB Antibody, FITC conjugated



PACO46328

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Receptor for ASNA1/TRC40-mediated insertion of tail-anchored (TA) proteins into the ER membrane.

**Gene ID:**

WRB

**Uniprot**

O00258

**Synonyms:**

Tail-anchored protein insertion receptor WRB (Congenital heart disease 5 protein) (Tryptophan-rich basic protein) (WRB), WRB, CHD5

**Immunogen:**

Recombinant Human Tail-anchored protein insertion receptor WRB protein (38-99AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A