

ZNF2 Antibody, Biotin conjugated



PACO46337

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in transcriptional regulation.

Gene ID:

ZNF2

Uniprot

Q9BSG1

Synonyms:

Zinc finger protein 2 (Zinc finger protein 2.2) (Zinc finger protein 661), ZNF2, ZNF661

Immunogen:

Recombinant Human Zinc finger protein 2 protein (134-289AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A