

atpB Antibody, Biotin conjugated



PACO46425

Product Information

Size:

50ug

Reactivity:

Arabidopsis thaliana

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Produces ATP from ADP in the presence of a proton gradient across the membrane. The catalytic sites are hosted primarily by the beta subunits.

Gene ID:

atpB

Uniprot

P19366

Synonyms:

ATP synthase subunit beta, chloroplastic (EC 3.6.3.14) (ATP synthase F1 sector subunit beta) (F-ATPase subunit beta), atpB

Immunogen:

Recombinant Arabidopsis thaliana ATP synthase subunit beta, chloroplastic protein (92-393AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A